

# Video Magnifier Solutions



Opti Verso



Quicklook 2



Prisma AF

 **ESCHENBACH OPTIK**  
Better vision. Better life.



[www.eschenbach.com](http://www.eschenbach.com)  
(800) 487-5389

# Video Magnifiers

■ Portable Hand-held

## ■ Portable Hand-helds

Quicklook  
Classic



### Quicklook Classic ASH-511-301

Convenient and compact, the Quicklook Classic is a cost-efficient, portable and lightweight video magnifier. Take it with you to read price tags, menus, train schedules, etc., and write checks or brief notes by flipping the camera into its writing position.

**Magnification:** 2x - 6x

**Viewing Modes:** Full color, enhanced b/w & w/b, with semi-colors

**Illumination:** LED

**Display:** 4.3" color TFT LCD

**Dimensions (l x w x h):** 6.7" x 3.7" x 1.1"

**Weight:** 11 oz.

**Power:** Built-in rechargeable Li-ion battery. 4-7 hour run time, 3 hour recharge time.

**Other:** Includes freeze frame feature. Camera rotates outward to allow for writing. Easy to grip with its ergonomic, rubberized casing.

Quicklook  
Zoom



### Quicklook Zoom ASH-511-201

The Quicklook Zoom is similar to the Quicklook Classic, but with a few additional features such as a larger magnification range and more viewing modes.

**Magnification:** 3x - 18x in 9 steps

**Viewing Modes:** Full color and high contrast w/ 56 color options and on-screen menu.

**Illumination:** LED

**Display:** 4.3" 16:9 high brightness LCD

**Dimensions (l x w x h):** 6.7" x 3.7" x 1.1"

**Weight:** 8.4 oz.

**Power:** Built-in rechargeable Li-ion battery. 4-7 hour run time, 3 hour recharge time.

**Other:** Includes freeze frame feature and on screen menu. Camera rotates outward to allow for writing. Easy to grip with its ergonomic, rubberized casing.

Quicklook  
Focus



### Quicklook Focus ASH-511-401

The Quicklook Focus is the premium product in the Quicklook line. In addition to all of the features of the Quicklook Zoom, the Quicklook Focus offers auto-focus and the ability to view objects up to 3 feet away.

**Magnification:** 3x - 18x in 9 steps

**Viewing Modes:** Full color and high contrast w/ 56 color options and on-screen menu.

**Illumination:** LED

**Display:** 4.3" 16:9 high brightness LCD

**Dimensions (l x w x h):** 6.7" x 3.7" x 1.1"

**Weight:** 8.4 oz.

**Power:** Built-in rechargeable Li-ion battery. 4-7 hour run time, 3 hour recharge time.

**Other:** Includes freeze frame feature and on screen menu. Camera rotates outward to allow for writing. Auto-focus feature allows for the viewing of objects a short distance away. Easy to grip with its ergonomic, rubberized casing.

## Portable Hand-helds

### Quicklook Touch ASH-511-601

The Quicklook Touch contains all of the attributes of the Quicklook Focus, but it also features new touch-screen technology. Up to 5 different touch zones are available on the screen giving users the ability to adjust magnification higher and lower, activate the auto-focus option, or access a menu to select other viewing options with a simple touch to the display screen.

**Magnification:** 3x - 18x in 9 steps

**Viewing Modes:** Full color and high contrast modes with 56 color options and an on-screen menu.

**Illumination:** LED

**Display:** 4.3" 16:9 high brightness LCD

**Dimensions (l x w x h):** 6.7" x 3.7" x 1.1"

**Weight:** 8.4 oz.

**Power:** Built-in rechargeable Li-ion battery. 4-7 hours run time, 3 hours recharge time.

**Other:** Includes freeze frame feature and on screen menu. Camera rotates outward to allow for writing notes. Auto-focus feature allows for viewing objects a short distance away. Easy to grip with its ergonomic, rubberized casing.



Quicklook Touch

### Quicklook 2 ASH-511-501

The Quicklook 2 has two built-in cameras integrated into one modern, ergonomic housing. Whether reading with the center mounted camera, or writing using the flip-out camera, or in hand-held mode, the Quicklook 2 enables you to find text and images easily and quickly. Perfect for reading menus, bills, schedules, phone numbers, ATM displays, price tags, food ingredients, prescription labels, maps or doing crossword puzzles, filling in forms.

**Magnification:** 3x - 24x when flat on text surface; 1x-6x when hand-held

**Viewing Modes:** Full color, positive, negative; select from 56 color combinations.

**Illumination:** LED

**Display:** 4.3" LCD. High brightness and high contrast ratio

**Dimensions (l x w x h):** 5.5" x 3.2" x .9"

**Weight:** 7.7 oz.

**Power:** Built-in rechargeable Li-ion battery with approximately 5-hours continuous use. 2.5 hours recharge time.

**Other:** Reading & writing are easily done with 2 integrated cameras.



Quicklook 2

### TVi Color ASH-508-103

The TVi Color is a hand-held, electronic magnifier that is plugged into your TV to magnify reading material. The camera is housed in the mouse-like device and can be used anywhere there is electricity. It features a handwriting stand for writing notes or signing checks.

**Magnification:** 9x - 18x (on a 14" TV), 13x - 26x (on a 21" TV)

**Viewing Modes:** full color, enhanced b/w, w/b

**Illumination:** LED

**Dimensions (l x w x h):** Camera is 5" x 2" x 2.2" in., base unit is 4" x 3" x 1.4" in.

**Weight:** Camera alone, 3.8 oz (total with base unit is 1 lb. 1 oz.)

**Power:** Electric

**Other:** The camera can be put in 'overview mode' for easy page orientation or can be turned around for use as a 'magnified mirror' for shaving or putting on make-up.



TVi Color

# Video Magnifiers

## ■ Portable Desktops

## ■ Portable Desktops

1654-104



### 1654-104 Color ELG Video Magnifier

With only 3 buttons, Eschenbach's Color ELG doesn't require reading a training manual before using it. The ELG connects to a TV and provides between 5x - 58.5x magnification, depending on the size of the TV screen. It provides 3 different viewing options: color, black and white, and reverse polarity (white on black). The X-Y table is padded to prevent reading materials from sliding, and the open system allows the user to view all sizes of magazines, books, prescription bottles, etc.

**Magnification:** 5x - 58.5x depending on TV size (TV not included)

**Viewing Modes:** Color, black on white and white on black

**Illumination:** Clip on fluorescent bulb included

**Dimensions (l x w x h):** 17" x 17" x 16.5"

**Weight:** 17 lbs.

**Power:** Electric

**Other:** Sliding X-Y Table Size: 17" x 13"

myReader2



### myReader2 by Humanware

The myReader2 is the only video magnifier that can be set by the user to scroll through text manually or automatically at varying speeds that suit the reader. Text can be arranged in column layout, row, word, or original layout. With its built-in 15" display, users can capture up to 13 pages at one time for reading continuity. There are 8 different viewing modes and the myReader2 folds up into a compact case for easy transportability.

**Magnification:** 0.7x - 45x

**Viewing Modes:** Color, black on white, white on black and 56 possible color combinations.

**Illumination:** 2 fluorescent and 1 halogen bulbs

**Display:** 15" TFT LCD display

**Dimensions (l x w x h):** 19" x 14.5" x 20"

**Weight:** 22 lbs.

**Power:** Electric

**Other:** Includes footswitch, remote control panel,

Fusion



### Fusion ASH-515-001

The Fusion is a portable video magnifier with a built-in monitor that rolls over text and images and magnifies them. With its tiltable display and easily removeable mouse-style camera, it is an excellent solution for users that need magnification for near.

**Magnification:** 6x - 16x with high brightness display

**Viewing Modes:** Full color, enhanced positive and negative reading modes

**Illumination:** LED

**Display:** 7 in. TFT LCD display

**Dimensions (l x w x h):** 7" x 7" x 2"

**Weight:** 2 lb. 4 oz.

**Power:** Built-in rechargeable Li-ion battery. 3 hour run time, 4 hour recharge time.

**Other:** Fusion can be used while it is recharging. Legs on camera allow it to be sit on an angle for writing. Carrying case included.

### Portable Desktops

#### Prisma ASH-505-002

The Prisma is a simply designed video magnifier that connects to a TV to view text or images. The Prisma folds into a compact, highly transportable configuration weighing just 2 lbs. 10 oz. Perfect for hobbies, such as needlework, making models, drawing and painting. The Prisma features brightness controls, LED illumination, a writing mat to prevent slippage, and carrying case.

**Magnification:** 3.8x - 35x (on a 14" TV screen)

**Viewing Modes:** Full color display mode with black on white and white on black modes.

**Illumination:** LED

**Dimensions (l x w):** 16" x 9"

**Weight:** 2 lbs. 10 oz.

**Power:** Electric

**Other:** Includes a carrying case and a writing mat that facilitates ease of writing with notes, crossword puzzles, bills, and letters.



Prisma

#### Prisma AF ASH-521-002

Prisma AF introduces a new level of simplicity and ease of use with the integration of full autofocus technology. Enjoy the broad range of magnification levels and premium image quality with the touch of a button. Premium image quality decreases eye fatigue; users can read more quickly and comfortably for longer periods of time.

**Magnification:** 2.5x - 30x (on a 17" TV screen). Greater magnification is available on larger screens.

**Viewing Modes:** Full color display mode with black on white and white on black modes.

**Illumination:** LED - the direction can be varied

**Dimensions (l x w x h):** 5.5" x 3.1" x .9" when flat

**Weight:** 2 lbs. 8 oz.

**Power:** Electric

**Other:** Includes a carrying box and a writing mat that facilitates ease of writing with notes, crossword puzzles, bills, and letters.



Prisma AF

#### Liberty Solo ASH-509-401

The Liberty Solo is an ergonomically designed, battery powered, portable video magnifier with a 12" screen. Choose a comfortable reading position and enjoy its tiltable display and lightweight mouse camera. Writing is also possible with the flip-out stand on the mouse camera.

**Magnification:** 1x - 30x

**Viewing Modes:** Full color, enhanced positive and negative reading modes, black on white and white on black and a range of semi colors, which are controlled by ergonomic buttons on the mouse unit.

**Illumination:** LED

**Display:** 12" TFT LCD

**Dimensions (l x w x h):** 12" x 10" x 2.4"

**Weight:** 4 lb. 10 oz. (mouse camera: 3.5 oz.)

**Power:** Built-in rechargeable Li-ion battery. 2 hour run time, 3 - 4 hour recharge time.

**Other:** Includes line markers that facilitate reading. Mouse camera features a handwriting stand that flips out to facilitate ease of writing. Carrying case and charger included.



Liberty Solo

#### Opti Lite ASH-516-103

Combine the freedom of your laptop with the Opti Lite to create a versatile, mini desktop magnification solution. Ideal for reading and writing, the Opti Lite has a wide range of magnification, fine focus control, and a 10-hour battery life. It's lightweight (2lbs., 10 oz.), so it can be easily carried from one location to another.

**Magnification:** 4x - 40x

**Viewing Modes:** Full color, complete range of enhanced color options with 10 customizable views.

**Illumination:** LED

**Dimensions (l x w):** 16" x 9"

**Weight:** 2 lbs. 10 oz.

**Power:** Built in rechargeable Li-ion battery.

**Other:** Third generation of Opti software is included and provides unique toolbar interface, fully adjustable image window size, image capture and save facility, and full control of video display attributes. Protective carrying case included.



Opti Lite

# Video Magnifiers

## Transportable Desktops/Desktops

Presto



Eclipse



iFlex



Horizon



## Transportable Desktops

### Presto ASH-517-001

The Presto is our transportable desktop video magnifier. Its plug-and-play design makes it easy to use; by simply connecting the power cord, it's ready to go. No assembly or set up is required!

**Magnification:** 2.5x - 40x

**Viewing Modes:** Color, black on white and white on black.

**Illumination:** LED

**Display:** 15" TFT LCD

**Dimensions (l x w x h):** 15" x 18" x 7" (when monitor is pulled up, height expands to 21" high)

**Weight:** 20 lbs. 15 oz.

**Power:** Electric

**Other:** Large XY table with breaking system accommodates full page viewing. Features include auto-focus camera and buttons to zoom-in and zoom-out, change image mode and change image enhancement. Optional storage case available with wheels and telescoping handle.

### Eclipse ASH-512-002

The Eclipse is similar to the Presto, but it includes additional features such as enhanced viewing modes, line markers and blinds. The Eclipse is a powerful performer at home, work, or in an educational environment.

**Magnification:** 1.5x - 32x

**Viewing Modes:** Color, black on white and white on black and a wide range of semi colors.

**Illumination:** LED

**Display:** 15" TFT screen, ergonomically positioned for comfort.

**Dimensions (l x w x h):** 15" x 18" x 7" (when monitor is pulled up, height expands to 21" high)

**Weight:** 20 lbs. 15 oz.

**Power:** Electric

**Other:** Large XY table with breaking system accommodates full page viewing. Features include auto-focus camera and buttons to zoom-in and zoom-out, change image mode and change image enhancement. Optional storage case available with wheels and telescoping handle.

### iFlex ASH-506-501

The iFlex is a unique solution for magnification at home or work where maximizing desk space is key. Easy to attach and detach from a standard computer monitor arm (not included), it allows for easy transportation.

**Magnification:** 3x - 50x

**Viewing Modes:** Color, black on white and white on black and a wide range of semi colors.

**Illumination:** LED

**Display:** 15" TFT LCD

**Dimensions (l x w x h):** 17.5" x 16.5" x 4.5"

**Weight:** 9 lbs.

**Power:** Electric

**Other:** Customized carrying case, footswitch, and video lead included.

## Desktops

### Horizon ASH-518-002

The Horizon is a new video magnifier with a host of built-in features, such as a large 19" wide screen (with anti-glare coating), tiltable monitor (so users can read comfortably for extended periods of time), auto and manual focus, and adjustable line markers and shutters that show users exactly where they are on the page at all times.

**Magnification:** 3.5x - 50x

**Viewing Modes:** Full color and a wide range of semi-colors

**Illumination:** LED

**Display:** 19" TFT LCD

**Dimensions (l x w x h):** 16.5" x 18" x 21"

**Weight:** 23.5 lbs.

**Power:** Electric

**Other:** Large XY table with breaking system accommodates full page viewing. Features both automatic and manual focus.

## Portable Desktops w/ Distance Camera

### Liberty Scholar ASH-509-202

An ideal magnification solution that includes a mouse camera, 12" integrated display, and a top-mounted pan-and-tilt zoom camera. Rechargeable batteries allow the system to be used when there is no access to power outlets.

**Magnification:** 12x - 27x

**Viewing Modes:** Full color and enhanced display modes

**Illumination:** LED

**Display:** 12" TFT LCD

**Dimensions (l x w x h):** 11.8" x 9.8" x 2.4"

**Weight:** 5 lb. 2 oz.

**Power:** Built-in rechargeable Li-ion battery. 1-3 hour run time, 3 - 4 hour recharge time.

**Other:** Distance camera offers 10x optical zoom.

### Opti Verso ASH-516-203

The ultimate portable reader for professionals and students who require distance AND close-up viewing, as well as writing capability. The Opti Verso has a collapsible design for easy portability.

**Magnification:** 2x - 20x (80x digital) on a typical 15" screen

**Viewing Modes:** Full color, range of color options, 10 customizable views

**Dimensions (l x w x h):** 16" x 9.4" x 13" (collapsing to 2.5")

**Weight:** 2 lbs. 11 oz.

**Power:** Built-in rechargeable Li-ion battery. 6 hour run time, 3 - 4 hour recharge time.

**Other:** USB lead, Opti software, power cord, power supply/charger, and carrying case included. 10x zoom control is possible from the camera or keyboard.

## Transportable Desktops w/ Distance Camera

### i-Qon ASH-519-101

The i-Qon is a 2-piece system that allows the TVi Color or Prisma to be used with a computer monitor. Split screen technology lets users view both their computer screen and the image projected by the TVi Color or Prisma on their computer monitor. The i-Qon offers adjustable line markers which facilitate easier text tracking and reading.

**Camera Image Modes:** Full color with complete range of enhanced color options.

**Weight:** Main unit - 17 oz.; keypad - 11 oz.

**Dimensions (l x w x h):** Main unit - 8.8" x 4" x 1"; Key pad - 4" diameter, 1" in height.

**Other:** Main unit, keypad, power supply and cable, video lead, and footswitch included. (A camera is required for use with this system.)

### Agilo ASH-520-002

The Agilo is a 3-piece system (the 2-piece i-Qon and a camera mounted on an adjustable arm) that works with a computer and existing monitor. Split screen technology lets users simultaneously view both their computer screen and the image projected by the camera on their computer monitor.

**Camera Image Modes:** Full color and high-contrast reading modes with a complete range of color options. Features line and line marker modes, and freeze image function.

**Weight:** Main unit - 17 oz.; keypad - 11 oz.

**Other:** Main control unit, keypad, camera w/ adjustable arm, power supply and cable, video lead, and foot switch included.

### Eclipse Scholar ASH-512-202

The Eclipse Scholar is the premier one-piece integrated design that can be used for reading, writing, and viewing distance objects. For use in the home, at work, or in an educational environment, it contains a top-mounted camera, full color, continuous autofocus, line markers and overview mode.

**Magnification:** 1.5x - 30x, 10x zoom on distance camera

**Viewing Modes:** Color, black on white and white on black and a wide range of semi colors.

**Illumination:** LED

**Display:** 15" TFT screen, ergonomically positioned for comfort

**Dimensions (l x w x h):** 15" x 18" x 7" (when monitor is pulled up and camera is attached, height expands to 25" high)


**Weight:** 21 lbs. 7 oz.

**Power:** Electric


**Other:** Large XY table with breaking system accommodates full page viewing. Features auto-focus camera with mirror and freeze modes. Camera, mounting bracket, Eclipse, carrying case w/wheels, and all power supplies included.

## Video Magnifiers

### Desktops w/ Distance Camera



Liberty  
Scholar



Opti  
Verso



i-Qon



Agilo



Eclipse  
Scholar

# Video Magnifier Solutions



Fusion



Agilo



Liberty Scholar

 **ESCHENBACH OPTIK**  
Better vision. Better life.

Eschenbach Optik of America, Inc.  
904 Ethan Allen Highway  
Ridgefield, CT 06877  
Phone (800) 487-5389  
Fax (888) 799-7200  
www.eschenbach.com  
Email: info@eschenbach.com

- ©2010 - All Rights Reserved.
- ESCHENBACH® is a Registered Trade Mark of ESCHENBACH OPTIK GmbH + Co

