



In this newsletter:

- QuickLook 2 – The Revolutionary Combination of Design and Technology!
- ‘SportScope’ Promotion – Buy One, Get One Free! Details Below
- Meet Your New Territory Manager – If you practice in MN, IA, IL or WI
- Free Low Vision Webinars . . . New Webinars Have Been Added – Details Below
- Eschenbach Will Be Exhibiting at These Tradeshows

QuickLook 2



*Revolutionary hand held magnification! The Perfect Combination of Design and Technology!
Dual cameras and touch-screen capability integrated for optimum ease of use whether reading text
that is near or at a distance, or for writing.*

The unique Quicklook 2 has two built-in cameras integrated into one modern, ergonomic housing. Whether reading with the center-mounted camera, writing under the flip out camera, or using the device in hand-held mode, the Quicklook 2 will enable you to find text and images quickly and easily.

With both manual and electronic focus, a superb depth of field is achieved. The Quicklook 2 also allows you to control and arrange images/text and select magnification that suits your patient’s specific needs.

This product comes in a sturdy case and has a robust design that will make it last for many years while delivering a ‘vision for living’.



Features and Specifications:

- SCREEN:** 4.3” LCD illumination, producing high brightness and contrast ratio
- MAGNIFICATION RANGE:** 3X to 24X (when flat on text surface) <1X to 6X (when handheld)
- MODES:** Full-color, Positive, Negative modes. User can select from 56 color combinations
- WEIGHT:** 220g or 7.7 ounces

UNIT SIZE: 140 x 82 x 23mm or 5.5” x 3.2” x .9”

POWER: Li-ion battery with approximately 5 hours continuous use. Recharge time is 2.5 hours.

READING / WRITING: Integrated reading and handwriting camera

To purchase this product, contact your Territory Manager or our Customer Service Department at (800) 487-5389.

Looking for a 'SportScope' for your Patients? Eschenbach offers the best 3x binocular mounted telescope for distance viewing with this special offer:


Buy One – Get One FREE!

Buy Two – Get Two Free, Buy Three – Get Three FREE, etc.

The 1634 Galilean Binocular System for Distance series offers you and your patients the most adjustable, comfortable, durable optics available, as well as the superior Eschenbach optics you've come to expect!

1634 Binocular System Features include:

- Magnification: 3x
- Visual Field: 9.5°
- PD range: 54 to 74 mm; easily adjustable and can be 'locked in'
- Lens: Achromatic objective, coated, 23 mm diameter
- Frame: Black plastic, reinforced with fiberglass
- Weight: 70 grams (2.4 ounces)
- Cold bend wire core temples
- Soft, metal reinforced saddle bridge (height and width adjustable)
- Telescope can be prescription-corrected
- Case: Simulated leather zippered case included



Take advantage of this perfect opportunity to build your inventory and begin offering a superior solution for your patient's distance viewing needs with a system that can be fitted immediately with no lab work required!

(This promotion runs thru June 30, 2010. Products purchased under this promotion may not be returned for credit at a later date.)

Minnesota, Iowa, Illinois and Wisconsin Meet Your New Territory Manager



Eschenbach is pleased to announce the following addition to our sales team: **Bill Stahl, Territory Manager for the North Central States**, has many years experience in the optical industry and will be calling on you soon. **To speak with Bill, call (800) 487-5389 x309. To contact Bill via email: bstahl@eschenbach.com**

We urge you to utilize your Territory Managers' practice management and product expertise! Your Territory Manager can help you grow your low vision care practice by providing in-services, consultative support, and by keeping you up-to-date on the newest technology and products that can benefit your low vision patients.

Free Low Vision Webinars!

Our webinars are 1-hour valuable learning experiences that cover a wide range of Low Vision topics that can be viewed and listened to from your office!

- Webinars are fast, convenient, and easy to attend! They are ideal for multiple listeners – you can use the speakerphone on your phone or your computer speakers so many of your staff can listen in at once.
- All webinars begin at 12pm Eastern time (9am Pacific) and last for 1-hour. At the end of each webinar, there will be time for questions.

Upcoming 2010 Webinars:

Date	Topic	Speaker
Wed., June 2 nd 2010	Low Vision Care-What's It All About?	Joe Chung, Territory Manager Southwest
Wed. June 16 th , 2010	Low Vision Clinical Tips & Tricks	Thomas I. Porter, OD, Director of Low Vision Services, Saint Louis University Dept. of Ophthalmology
Wed., July 7 th 2010	Video Magnifiers – Everything You Need to Know About Them	Ron Dare, Territory Manager Mid-Atlantic States
Wed., August 4 th 2010	Training the Patient on Low Vision Devices + What Should be Done BEFORE the First Visit	Eugene Pfeifer, Territory Manager Northeast

To register for webinars or for additional details, email dsommers@eschenbach.com or call Doris Sommers at (800) 487-5389 x359.



May / June / July 2010 Trade Show Calendar

Eschenbach will be exhibiting at these tradeshows. Visit our booth to see our newest products!



DATE	TRADESHOW	LOCATION
May 24-25, 2010	Mississippi AER	Flowood, MS
May 26, 2010	CNIB – Toronto	Toronto, ON CANADA
May 27-30, 2010	BCAO-AAO Conference - Vision Rocks 2010	Banff Springs, AB CANADA
June 3-4, 2010	Alabama School f/t Blind Technology Symposium	Talledega, AL
June 15-18, 2010	TX Assistive Technology Network Statewide Conf.	San Antonio, TX
June 16-19, 2010	American Optometric Association Conference	Orlando, FL

For more information on any of the above tradeshows, please check the Internet or contact your local Eschenbach Territory Manager at (800) 487-5389.

(To learn more about the products listed in this email including pricing information, or to arrange an appointment with one of our representatives, please call us at (800) 487-5389 or visit our website at www.eschenbach.com. This email is sent on a regular basis to low vision practitioners to provide product updates, dispensing tips and news from Eschenbach. If you do not wish to receive future emails, please send an email to dsommers@eschenbach.com and ask that your name be removed from our email list.)



ESCHENBACH OPTIK
Better vision. Better life.

Toll-free (800) 487-5389
www.eschenbach.com